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|                                                                                                                                     |         |
|-------------------------------------------------------------------------------------------------------------------------------------|---------|
| NIBP Cuff Reusable, Adult, Thigh, Single M5306 Bladderless, Tube Length 20 cm                                                       | 13-082  |
| NIBP Cuff Disposable, Neonate, Single M5541-1# with CT-167 Connector                                                                | 13-085  |
| NIBP Cuff Disposable, Neonate, Single M5541-2# with CT-167 Connector                                                                | 13-086  |
| NIBP Cuff Disposable, Neonate, Single M5541-3# with CT-167 Connector                                                                | 13-087  |
| NIBP Cuff Disposable, Neonate, Single M5541-4# with CT-167 Connector                                                                | 13-088  |
| PU Legthing Tube (Black)                                                                                                            | 13-097  |
| NIBP Cuff Reusable – Adult – Single<br>M5114PU, TPU Bladder, Tube Length 20 cm                                                      | 13083   |
| NIBP Cuff Reusable – Adult – Single<br>M5104 Nylon, TPU Bladder, Tube Length 20 cm                                                  | 13084   |
| <b>TEMP</b>                                                                                                                         |         |
| TEMP Probe, Skin ,LAUNCH (LNHmed) (98ME04GA634)                                                                                     | 10-083  |
| TEMP Probe,Rectal,LAUNCH (LNHmed) (98ME04GA635)                                                                                     | 10-084  |
| TEMP Interface Probe, Data Cable for Redel Connector to Temp Probe                                                                  | 24-073  |
| <b>IBP</b>                                                                                                                          |         |
| IBP Transducer, MEDEX , MX860/866 Novatrans                                                                                         | 16-001  |
| IBP Disposable Dome , MEDEX , MX860/866 Novatrans Dome                                                                              | 16-031  |
| IBP Extension Cable , MEDEX , MX860/866 Novatrans Extension                                                                         | 16-042  |
| IBP Transducer, MEDEX , MX960 Logical                                                                                               | 16-002  |
| IBP Disposable Dome, MEDEX , MX960 Logical Dome                                                                                     | 16-033  |
| IBP Extension Cable , MEDEX , MX960 Logical Extension                                                                               | 16-043  |
| IBP Transducer Cable , TRUWAVE                                                                                                      | 16-037  |
| IBP Transducer , Disposable ,RX only ,PX260                                                                                         | 16-036  |
| IBP Interface Probe, One channel IBP interface                                                                                      | 16-051  |
| IBP Interface Probe, Two channel IBP interface                                                                                      | 16-052  |
| IBP Transducer kit, Disposable, iPex, Ref BKT,164ET                                                                                 | 16-046  |
| IBP Cable, Ipex, P/N: BKT,164ET                                                                                                     | 16-053  |
| IBP Bracket for iPex Transducer                                                                                                     | 16-047  |
| <b>CO2 (main stream)</b>                                                                                                            |         |
| IRMA Disposable Airway Adapter without O2 port                                                                                      | 20-025  |
| IRMA Disposable Airway Adapter for infant                                                                                           | 20-035  |
| IRMA Adapter Cable                                                                                                                  | 24-111  |
| Probe Holder for IRMA sensor                                                                                                        | 20-043  |
| CO2Airway daptor,Disposable,neonate/pediatric                                                                                       | 20091   |
| CO2Airway Adaptor,Disposable,Adult/pediatric                                                                                        | 20092   |
| <b>CO2 (side stream)</b>                                                                                                            |         |
| Nomoline with luer lock connector. 2 m. Box of 25                                                                                   | 20-045  |
| Clamp of ISA Module Holder                                                                                                          | 20-055  |
| VersaStream,CO2/GasAirway,Adapter Sampling Line, Adult / Pediatric                                                                  | 20-077  |
| VersaStream,CO2/GasAirway,Adapter Sampling Line, Infant                                                                             | 20-078  |
| VersaStream, CO2/Gas Sampling Line with Luer Lock Male(it uses with Sidestream Airway Adapter-Adult/Pediatric, part number:4420531) | 20-079  |
| T4F Water Filter for Capno-S+                                                                                                       | 20-094  |
| Sample line for Capno-S+                                                                                                            | 20-095  |
| T Airway Adapter for Capno-S+                                                                                                       | 20-096  |
| <b>BFA</b>                                                                                                                          |         |
| BFA Accessory Patient Cable, SAADAT                                                                                                 | 22- 028 |
| <b>Adapter</b>                                                                                                                      |         |
| Saadat Adaptor 60W, 15v for Aria                                                                                                    | 09263   |
|                                                                                                                                     |         |

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| ECG Electrodes                                                                                                                                                                                           | Part num.                     |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------|
| Adults ECG Disposable Electrodes, FIAB Manufacturer                                                                                                                                                      | P28042<br>(REF: F9060)        |
| Pediatric ECG Disposable Electrodes, FIAB Manufacturer                                                                                                                                                   | P28047<br>(REF: F9060P)       |
| Arbo H124SG, COVIDIEN Manufacturer                                                                                                                                                                       | P10079<br>(REF: 31.1245.21)   |
| EEG Electrodes                                                                                                                                                                                           |                               |
| Neuroline 720, AMBU Manufacturer                                                                                                                                                                         | P22009<br>(REF: Neuroline720) |
| Headset for Aria with TC Station                                                                                                                                                                         |                               |
| Headset,A4TECH,HS-5P for Aria+ TC Station<br>Headphone Specification:<br>●Frequency response:20Hz-20khz<br>●Impedance:32ohm<br>●Sensitivity:97Db/10mw/3mm,1mwJack pulg:OD3.5mm<br>●StereoCable Length:2m | P09322                        |

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| CLASSIFICATION                  |                                                                                                                   |
|---------------------------------|-------------------------------------------------------------------------------------------------------------------|
| Protection against electroshock | Class I, Type CF for all modules (except CO2 module & NIBP module that are BF) (based on IEC 60601-1)             |
| Mode of operation               | Continuous operation equipment                                                                                    |
| Harmful Liquid Proof Degree     | Aria monitor: IP32<br>Stations & Adaptor: IPX1                                                                    |
| Method of disinfection          | Refer to each module's chapters and chapter Care Cleaning for detail.                                             |
| Safety of anesthetic mixture    | Not suitable for use in the presence of a flammable anaesthetic mixture with air or with oxygen or nitrous oxide. |

| GENERAL            |                                                                                                                                                                                                                                                                              |
|--------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Display            | COLOR TFT 480 × 272" 5" Flexible display Configuration                                                                                                                                                                                                                       |
| Waveforms          | ECG, SPO2, RESP/CO2 (Aria-TC: Freezable), IBP1,IBP2,IBP3,IBP4, EEG(Aria: Freezable)                                                                                                                                                                                          |
| Numeric Parameters | Aria:<br>HR, PVCs,ST, SPO2 (%SPO2, PR, PI), Rainbow (SpMet, SpCO, SpHb, SpOC, PVI), NIBP (SYS, DIA, MAP), RR, TEMP1,2, IBP1,2,3,4 (SYS, DIA, MAP), EtCo2, FiCo2, AWRR, BFI, BS%, EMG%, SQI%,<br><br>Aria-TC:<br>HR, SpO2 (%SpO2, PR),ST,PVCs, NIBP (SYS, DIA, MAP), RR, TEMP |
| Operation Method   | Membrane, Touch screen                                                                                                                                                                                                                                                       |
| AC Power(Adaptor)  | 100 - 240 VAC, 50/60 Hz , Ip: 1.4 – 0.7 A, Output:15VDC,4A                                                                                                                                                                                                                   |
| Application        | Compact and Mobile Monitor.                                                                                                                                                                                                                                                  |
| Safety             | Based on IEC 60601-1, Class I                                                                                                                                                                                                                                                |
| Protection         | Against Electro surgery and Defibrillator and EMC.                                                                                                                                                                                                                           |

| ECG              |                                                             |
|------------------|-------------------------------------------------------------|
| Leads            | Selectable: 3 ,5 or 10 Wires                                |
|                  | For 3 wire: I, II, III                                      |
|                  | For 5 wire :I,II,III,V,aVR,aVF,aVL                          |
|                  | For 10 wire : I,II,III, aVR,aVF,aVL , V1,V2, V3, V4, V5. V6 |
| Dynamic Range    | ± 5 mV                                                      |
| Lead Off Current | < 90 nA                                                     |
| Gain             | 4, 2, 1, 1/2, 1/4, Auto                                     |
| Calibration      | 1mV, 0.5 sec                                                |
| Filters          | “MONITOR” ( 0.5 - 24 Hz )                                   |
|                  | “NORMAL” ( 0.5 - 40 Hz )                                    |
|                  | “EXTENDED” ( 0.05-100 Hz)                                   |

|                           |                                                                             |                                                   |
|---------------------------|-----------------------------------------------------------------------------|---------------------------------------------------|
| CMRR                      | > 98 dB                                                                     |                                                   |
| Internal Noise            | < 30 $\mu$ V RTI                                                            |                                                   |
| Input Impedance           | > 5 M $\Omega$                                                              |                                                   |
| QRS Detection             | Duration                                                                    | 40 to 120 msec                                    |
|                           |                                                                             | 0.25 to 5 mV for Adult/Pediatric                  |
|                           | Amplitude                                                                   | 0.2 to 5 mV for Neonate                           |
| Heart Rate Range          | 15 - 300 BPM for adult/Pediatric                                            |                                                   |
|                           | 15 - 350 BPM for neonate                                                    |                                                   |
| Accuracy                  | $\pm$ 1% or 2 BPM                                                           |                                                   |
| Tall T-Wave               | Reject up to 1.2 mV Amp.                                                    |                                                   |
| Pacer Detection/Rejection | Duration                                                                    | 0.1 - 2 msec                                      |
|                           | Amp                                                                         | $\pm$ 2 to $\pm$ 700 mV (Without over/undershoot) |
|                           | Reject from heart rate counter.                                             |                                                   |
|                           | Re-insert into ECG to display on screen.                                    |                                                   |
|                           | Ineffective pace rejection                                                  | HR:0, Pace: 60                                    |
|                           |                                                                             | HR:60, Pace:60                                    |
|                           |                                                                             | HR:30, Pace:80                                    |
|                           | Beside rejection of atrial paces precede ventricular paces by 150 or 250 ms |                                                   |
| Protection                | Defibrillator and Electrosurgery                                            |                                                   |

| ST ANALYSIS        |                                                                |
|--------------------|----------------------------------------------------------------|
| Display resolution | 0.01 mV                                                        |
| Measurement Range  | -2mv to +2mv                                                   |
| Alarm Range        | -2mv to +2mv                                                   |
| Features           | User Adjustable Isoelectric and ST point trending of ST values |
| Update period      | 5 Sec.                                                         |

| ARRHYTHMIA ANALYSIS |                                                                                                    |
|---------------------|----------------------------------------------------------------------------------------------------|
| Type                | ASYS, VFIB, VTAC, RUN, AIVR, COUPLET, BIGEMINY, TRIGEMINY, TACHY, BRADY, AFIB, PAUS, FREQUENT PVCs |
| Learning            | Rapid Learning: only 20 seconds required for recognition of dominant rhythm.                       |
| Method              | Real time arrhythmia detection with innovative feature.                                            |
| Memory              | Capability of storing the latest 150 ARR event (waveform and Parameters)                           |

| <b>NIBP</b>                      |                                                                                                            |                   |
|----------------------------------|------------------------------------------------------------------------------------------------------------|-------------------|
| Physical dimensions              | Width: 2.52 in (64 mm)<br>Length: 3.11 in (79 mm)<br>Height: 0.96 in (24.5 mm)<br>Weight: 0.25lbs (0.11kg) |                   |
| Operating voltage                | 5V DC $\pm$ 5% Input Voltage                                                                               |                   |
| Current Consumption              | Idle: 52mA<br>Measurement: 190mA<br>Inflation: 350 mA                                                      |                   |
| Communication Protocol           | Serial RS232 or TTL                                                                                        |                   |
| Operating Conditions             | Temperature: 50C to 400C<br>Humidity: 20% to 90%                                                           |                   |
| Storage Conditions               | Temperature: -250C to 600C<br>Humidity: 10% to 100%                                                        |                   |
| Technique                        | Oscillometry                                                                                               |                   |
| Measurement time                 | 20-25 sec (excluding cuff inflation time)                                                                  |                   |
| Safety & Regulatory Standards    | IEC 60601-1, AAMI SP10, EN1060-1, EN 1060-3, EN1060-4, IEC 80601-2-30, ISO 81060-2: 2014                   |                   |
| NIBP Accuracy                    | Meets BS EN ISO 81060-2:2014                                                                               |                   |
| Measurement method               | Oscillometric                                                                                              |                   |
| Mode of operation                | Manual/Automatic/Stat                                                                                      |                   |
| Auto mode repetition intervals   | 1,2,3,5,10,15,20, 30, 45, 60, 90 minutes and 2, 4, 8,12,16,20 and 24 hours.                                |                   |
| STAT mode cycle time             | 5 min                                                                                                      |                   |
| Max measurement time             | Adult, pediatric: 180 s<br>Neonate: 90 s                                                                   |                   |
| Heart rate range                 | 30 to 240 bpm                                                                                              |                   |
| Cuff pressure Rang               | Adult                                                                                                      | 0-290 mmHg        |
|                                  | Neonate                                                                                                    | 0-145 mmHg        |
| Measurement Range(mmHg)          | Adult                                                                                                      | SYS 30 ~ 255 mmHg |
|                                  |                                                                                                            | DIA 15 ~ 220 mmHg |
|                                  |                                                                                                            | MAP 20 ~ 235 mmHg |
|                                  | Neonate                                                                                                    | SYS 30 ~ 135 mmHg |
|                                  |                                                                                                            | DIA 15 ~ 110 mmHg |
|                                  |                                                                                                            | MAP 20 ~ 125 mmHg |
|                                  | Pediatric                                                                                                  | SYS 30 ~ 240mmHg  |
|                                  |                                                                                                            | DIA 15 ~ 220 mmHg |
|                                  |                                                                                                            | MAP 20 ~ 230 mmHg |
| Software overpressure protection | Adult: 290 $\pm$ 3 mmHg<br>Pediatric: 240 $\pm$ 3 mmHg<br>Neonate: 145 $\pm$ 3 mmHg                        |                   |
| Resolution                       | 1 mmHg                                                                                                     |                   |
| Transducer accuracy              | $\pm$ 3 mmHg over full range in operating conditions                                                       |                   |

|                                                |                                                            |
|------------------------------------------------|------------------------------------------------------------|
| default initial cuff inflation pressure (mmHg) | Adult : 150 mmHg<br>Pediatric: 140mmHg<br>Neonate: 85 mmHg |
| Memory                                         | 100 Records                                                |

### SpO2 (Masimo Rainbow Set)

|                 |                             |                         |                                                    |
|-----------------|-----------------------------|-------------------------|----------------------------------------------------|
| SpO2 Parameters | SpO2,PI,PR                  |                         |                                                    |
| Method SpO2     | 2 Wavelengths of light used |                         |                                                    |
|                 | SpOC                        |                         |                                                    |
|                 | SpCO                        |                         |                                                    |
|                 | SpMet                       |                         |                                                    |
|                 | SpHb                        |                         |                                                    |
|                 | PVI                         |                         |                                                    |
| Method Rainbow  | 7+Wavelengths of light used |                         |                                                    |
| Range           | SpO2                        | 0 – 100 %               |                                                    |
|                 | PR                          | 25 – 240 bpm            |                                                    |
|                 | PI                          | 0.02– 20.0 %            |                                                    |
|                 | SpMet                       | 0 – 99.9 %              |                                                    |
|                 | SpCO                        | 0 – 99 %                |                                                    |
|                 | SpHb                        | 0 – 25.0 g/dL           |                                                    |
|                 | SpOC                        | 0 – 35.0 ml/dL          |                                                    |
|                 | PVI                         | 0 – 100 %               |                                                    |
| Accuracy        | Oxygen Saturation           |                         |                                                    |
|                 | No motion conditions        | Adult/Pediatric         | ±2% (SpO2 70 ~ 100%)<br>Unspecified (SpO2 0 ~ 69%) |
|                 |                             | Neonate                 | ±3% (SpO2 70 ~ 100%)<br>Unspecified (SpO2 0 ~ 69%) |
|                 | Motion conditions           | Adult/Pediatric/Neonate | ±3% (SpO2 70 ~ 100%)<br>Unspecified (SpO2 0 ~ 69%) |
|                 |                             |                         |                                                    |

|            |                              |                         |                                                    |
|------------|------------------------------|-------------------------|----------------------------------------------------|
|            | Low perfusion conditions     | Adult/Pediatric/Neonate | ±2% (SpO2 70 ~ 100%)<br>Unspecified (SpO2 0 ~ 69%) |
|            | Pulse Rate                   |                         |                                                    |
|            | No motion conditions         | Adult/Pediatric/Neonate | ±3bpm (PR 25 ~ 240)                                |
|            | Motion conditions            | Adult/Pediatric/Neonate | ±5bpm (PR 25 ~ 240)                                |
|            | Low perfusion conditions     | Adult/Pediatric/Neonate | ±5bpm (PR 25 ~ 240)                                |
|            | Carboxyhemoglobin Saturation |                         |                                                    |
|            | Carboxyhemoglobin Saturation | Adult/Pediatric         | ±3% (1 - 40)                                       |
|            | Methemoglobin Saturation     |                         |                                                    |
|            | Methemoglobin Saturation     | Adult/Pediatric/Neonate | ±1% (1 - 15)                                       |
|            | Total Hemoglobin             |                         |                                                    |
|            | Total Hemoglobin             | Adult/Pediatric         | ±1g/dL (8 - 17) g/dL                               |
| Resolution | SpO2                         | 1%                      |                                                    |
|            | PI                           | 0.1%                    |                                                    |
|            | PR                           | 1 BPM                   |                                                    |
|            | SpCO                         | 1.0 %                   |                                                    |
|            | SpMet                        | 0.1 %                   |                                                    |
|            | SpHb                         | 0.1 g/dL                |                                                    |
|            | PVI                          | 0-100%                  |                                                    |
|            | SpOC                         | 0.1 ml/dL               |                                                    |

1. The Masimo SET technology with Masimo sensors has been validated for no motion accuracy in human blood studies on healthy adult male and female volunteers with light to dark skin pigmentation in induced hypoxia studies in the range of 70-100% SpO<sub>2</sub> against a laboratory CO-Oximeter and ECG monitor. This variation equals  $\pm 1$  standard deviation. Plus or minus one standard deviation encompasses 68% of the population.
2. The Masimo SET technology with Masimo sensors has been validated for motion accuracy in human blood studies on healthy adult male and female volunteers with light to dark skin pigmentation in induced hypoxia studies while performing rubbing and tapping motions, at 2 to 4 Hz at an amplitude of 1 to 2 cm and a non-repetitive motion between 1 to 5 Hz at an amplitude of 2 to 3 cm in induced hypoxia studies in the range of 70-100% SpO<sub>2</sub> against a laboratory CO-Oximeter and ECG monitor. This variation equals  $\pm 1$  standard deviation, which encompasses 68% of the population.
3. The Masimo SET technology has been validated for low perfusion accuracy in bench top testing against a Biotek Index 2™ simulator and Masimo's simulator with signal strengths of greater than 0.02% and transmission of greater than 5% for saturations ranging from 70 to 100%. This variation equals  $\pm 1$  standard deviation. Plus or minus one standard deviation encompasses 68% of the population.
4. The Masimo SET Technology with Masimo Neo sensors has been validated for neonatal motion accuracy in human blood studies on healthy adult male and female volunteers with light to dark skin pigmentation in induced hypoxia studies while performing rubbing and tapping motions, at 2 to 4 Hz at an amplitude of 1 to 2 cm and a non-repetitive motion between 1 to 5 Hz at an amplitude of 2 to 3 cm in induced hypoxia studies in the range of 70-100% SpO<sub>2</sub> against a laboratory CO-Oximeter and ECG monitor. This variation equals  $\pm 1$  standard deviation. Plus or minus one standard deviation encompasses 68% of the population. 1% has been added to the results to account for the effects of fetal hemoglobin present in neonates.
5. The Masimo SET technology with Masimo sensors has been validated for pulse rate accuracy for the range of 25 -240 bpm in bench top testing against a Biotek Index 2™ simulator. This variation equals  $\pm 1$  standard deviation. Plus or minus one standard deviation encompasses 68% of the population.
6. See sensor directions for use (DFU) for complete application information. Unless otherwise indicated, reposition reusable sensors at least every 4 hours and adhesive sensors at least every 8 hours.
7. Sensor accuracy specified when used with Masimo technology using a Masimo patient cable for LNOP sensors, RD SET sensors, the LNCS sensors, or the M-LNCS sensors. Numbers represent A (RMS error compared to the reference). Because pulse oximeter measurements are statistically distributed, only about two-thirds of the measurements can be expected to fall within a range of  $\pm A_{rms}$  compared to the reference value. Unless otherwise noted, SpO<sub>2</sub> accuracy is specified from 70% to 100%. Pulse Rate accuracy is specified from 25 to 240 bpm.
8. Masimo M-LNCS, LNOP, RD SET, and LNCS sensors types have the same optical and electrical properties and may differ only in application type (adhesive/non-adhesive/hook & loop), cable lengths, optical component locations (top or bottom of sensor as aligned with cable), adhesive material type/size, and connector type (LNOP 8 pin modular plug, RD 15 pin modular plug, LNCS 9 pin, cable based, and M-LNCS 15 pin, cable based). All sensor accuracy information and sensor application instructions are provided with the associated sensor directions for use.

## TEMPERATURE

|            |                                                                          |
|------------|--------------------------------------------------------------------------|
| Channel    | Aria:<br>2 Channel<br>Monitoring 2 channel.<br><br>Aria-TC:<br>1 channel |
| Probe Type | YSI 400 Compatible                                                       |
| Range      | 0 - 50 °C                                                                |
| Accuracy   | $\pm 0.2$ °C                                                             |

## RESPIRATION

|                   |               |
|-------------------|---------------|
| Method            | Impedance     |
| Base Resistance   | 250 -1250 Ohm |
| Dynamic Range     | 0.2 - 2 Ohm   |
| Breath Rate Range | 0 - 253 BrPM  |

|          |               |
|----------|---------------|
| Accuracy | ±2% or 2 BrPM |
|----------|---------------|

| IBP                      |                                    |
|--------------------------|------------------------------------|
| Channel                  | 4 Channels                         |
| Measurement Range        | SYS -50 ~ 300 mmHg                 |
|                          | DIA -50 ~ 300 mmHg                 |
|                          | MAP -50 ~ 300 mmHg                 |
| Pressure Filter          | 8Hz, 16Hz, 22Hz selectable         |
| Press Sensor Sensitivity | 5 $\mu$ V / V / mmHg               |
| Resolution               | 1 mmHg                             |
| Accuracy                 | 2 % or 2mmHg (each one is greater) |

| Multi-gas, Mainstream (MASIMO SWEDEN AB)                                    |                                                                                                 |                                                |
|-----------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------|------------------------------------------------|
| IRMA CO2                                                                    | CO2                                                                                             |                                                |
| IRMA AX+                                                                    | CO2, N2O, primary and secondary agents (HAL, ISO, ENF, SEV, DES)                                |                                                |
| Gas /CO2 Interface                                                          | Connector and S/W Interface Driver, Applicable for All Gas and CO2 Modules.                     |                                                |
| Description                                                                 | Extremely compact infrared mainstream multigas probe available in two parameter configurations. |                                                |
| Cable length                                                                | 2.5 m $\pm$ 0.1 m                                                                               |                                                |
| Recovery time after defibrillator test                                      | Unaffected                                                                                      |                                                |
| Drift of measurement accuracy                                               | No drift                                                                                        |                                                |
| Surface temperature<br>(at ambient temp. 23°C)                              | IRMA CO2                                                                                        | Max 39°C / 102°F                               |
|                                                                             | IRMA AX+                                                                                        | Max 46°C / 115°F                               |
| Interface                                                                   | Modified RS-232 serial interface operating at 9600 bps.                                         |                                                |
| Airway adapters                                                             | Disposable adult/pediatric:                                                                     | - Adds less than 6 ml deadspace.               |
|                                                                             |                                                                                                 | - Pressure drop less than 0.3 cm H2O @ 30 LPM. |
|                                                                             | Disposable infant:                                                                              | - Adds less than 1 ml deadspace.               |
|                                                                             |                                                                                                 | - Pressure drop less than 1.3 cm H2O @ 10 LPM. |
|                                                                             | (Infant Airway Adapter recommended for Tracheal Tube ID size = 4 mm)                            |                                                |
| Degree of protection against harmful ingress of water or particulate matter | IPX4                                                                                            |                                                |

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |                                                                                                                                |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------|
| Method of sterilization                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      | The IRMA system contains no sterile parts.                                                                                     |
| Mode of operation                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            | CONTINUOUS OPERATION                                                                                                           |
| Data output                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |                                                                                                                                |
| Breath detection                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             | Adaptive threshold, minimum 1 vol% change in CO <sub>2</sub> concentration.                                                    |
| Respiration rate <sup>1</sup>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                | 0–150 ±1 bpm. The respiration rate is displayed after three breaths and the average value is updated every breath.             |
| <p>Fi and ET<sup>2</sup></p> <p>Fi and ET are displayed after one breath and have a continually updated breath average.</p> <p>The following methods are used to calculate end-tidal (ET) values:</p> <p>-CO<sub>2</sub>: The highest concentration of CO<sub>2</sub> during one breathing cycle with a weight function applied to favor values closer to the end of the cycle.</p> <p>-N<sub>2</sub>O and anesthetic agents: The momentary gas concentration at the time point where ET-CO<sub>2</sub> is detected.</p> <p>ET-values for anaesthetic agents and N<sub>2</sub>O (IRMA AX+) will typically decrease below nominal value when respiration rate exceeds 80 bpm. The maximum decrease is described by the formula <math>ET = 80 \cdot ET_{nom} / RR</math>.</p> <p>ET-CO<sub>2</sub> will be within specification for all respiration rates up to 150 bpm (IRMA AX+ and IRMA CO<sub>2</sub>)</p> |                                                                                                                                |
| Automatic agent identification                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               | IRMA AX+: Primary and secondary agent.                                                                                         |
| Gas Analyzer                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |                                                                                                                                |
| Probe                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        | 2-9 channel NDIR type gas analyzer measuring at 4-10 µm. Pressure, temperature and full spectral interference correction.      |
| Calibration                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  | Zeroing recommended when changing Airway adapter (IRMA AX+)<br><br>No span calibration required for the IR bench.              |
| Warm-up time                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 | IRMA CO <sub>2</sub> : < 10 seconds (concentrations reported and full accuracy)                                                |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              | IRMA AX+: < 20 seconds (concentrations reported, automatic agent identification enabled and full accuracy)                     |
| Rise time <sup>3</sup> (@10 l/min)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           | CO <sub>2</sub> ≤ 90 ms                                                                                                        |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              | N <sub>2</sub> O ≤ 300 ms                                                                                                      |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              | HAL, ISO, ENF, SEV, DES ≤ 300ms                                                                                                |
| Primary agent threshold                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      | 0.15 vol%. When an agent is identified, concentrations will be reported even below 0.15 vol% as long as apnea is not detected. |

Measured at I/E ratio 1:1 using breath simulator according to EN ISO 80601-2-55 fig. 201.101.<sup>1</sup>

Measured at I/E ratio 1:1 using breath simulator according to EN ISO 80601-2-55 fig. 201.101.<sup>2</sup>

Measured @ 10 l/min with gas concentration steps corresponding to 30% of total measuring range for each gas.<sup>3</sup>

|                                                                                                                                                                                                                                                        |                             |                                            |                             |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------|--------------------------------------------|-----------------------------|
| Secondary agent threshold                                                                                                                                                                                                                              |                             | 0.2 vol% +10% of total agent concentration |                             |
| Agent identification time                                                                                                                                                                                                                              |                             | <20 seconds. (Typically < 10 seconds)      |                             |
| Total system response time <sup>1</sup>                                                                                                                                                                                                                |                             | < 1 second                                 |                             |
| Accuracy - standard conditions                                                                                                                                                                                                                         |                             |                                            |                             |
| The following accuracy specifications are valid for dry single gases at 22 ± 5 °C and 1013 ± 40 hPa                                                                                                                                                    |                             |                                            |                             |
| Gas                                                                                                                                                                                                                                                    | Range                       |                                            | Accuracy                    |
| CO2                                                                                                                                                                                                                                                    | 0 to 15 vol%                |                                            | ±(0.2 vol% +2% of reading)  |
| N2O                                                                                                                                                                                                                                                    | 0 to 100 vol%               |                                            | ±(2 vol% +2% of reading)    |
| HAL,ISO,ENF                                                                                                                                                                                                                                            | 0 to 8 vol%                 |                                            | ±(0.15 vol% +5% of reading) |
| SEV                                                                                                                                                                                                                                                    | 0 to 10 vol%                |                                            | ±(0.15 vol% +5% of reading) |
| DES                                                                                                                                                                                                                                                    | 0 to 22 vol%                |                                            | ±(0.15 vol% +5% of reading) |
| Accuracy - all condition                                                                                                                                                                                                                               |                             |                                            |                             |
| The following accuracy specifications are valid for all specified environmental conditions except for interference specified in the table “Interfering gas effects” and the section “Effects from water vapor partial pressure on gas readings” below. |                             |                                            |                             |
| Gas                                                                                                                                                                                                                                                    | Accuracy                    |                                            |                             |
| CO2                                                                                                                                                                                                                                                    | ±(0.3 kPa + 4% of reading)  |                                            |                             |
| N2O                                                                                                                                                                                                                                                    | ±(2 kPa + 5% of reading)    |                                            |                             |
| Agents <sup>2</sup>                                                                                                                                                                                                                                    | ±(0.2 kPa + 10% of reading) |                                            |                             |
| Gas concentration conversion                                                                                                                                                                                                                           |                             |                                            |                             |
| Gas concentration is reported in units of volume percent. The concentration is defined as:                                                                                                                                                             |                             |                                            |                             |
| $\% \text{ gas} = \frac{(\text{Partial pressure of gas component})}{(\text{Total pressure of gas mixture})} * 100$                                                                                                                                     |                             |                                            |                             |
| The total pressure of the gas mixture is estimated by measuring the actual atmospheric pressure in the IRMA probe.                                                                                                                                     |                             |                                            |                             |
| Effects from water vapor partial pressure on gas readings                                                                                                                                                                                              |                             |                                            |                             |

<sup>1</sup> Measured according to EN ISO 80601-2-55.

<sup>2</sup> The accuracy specification for IRMA AX+ is not valid if more than two agents are present in the gas mixture. If more than two agents are present, an alarm will be set.

The effects of water vapor are illustrated by the examples in the following table. The two columns to the right show the relative error in displayed concentrations when adding or removing water vapor from the gas mixture, and referencing the measurement to dry gas conditions at actual temperature and pressure (ATPD) or saturated conditions at body temperature (BTPS).

| Temp [C] | RH [%] | P [hPa] | H2O<br>part.pres.<br>[hPa] | errrel [%] | errrel ATPD<br>[%] | errrel [%]<br>BTPS |
|----------|--------|---------|----------------------------|------------|--------------------|--------------------|
| 10       | 20     | 1013    | 2                          | 0          | -0.2               | +6.0               |
| 20       | 20     | 1013    | 5                          | 0          | -0.5               | +5.7               |
| 25       | 0      | 1013    | 0 (ATPD)                   | 0          | 0                  | +6.2               |
| 25       | 23     | 1013    | 7.3                        | 0          | -0.7               | +5.5               |
| 25       | 50     | 1013    | 16                         | 0          | -1.6               | +4.6               |
| 30       | 80     | 1013    | 42                         | 0          | -4.1               | +2.0               |
| 37       | 100    | 1013    | 63 (BTPS)                  | 0          | -6.2               | 0                  |
| 37       | 100    | 700     | 63                         | 0          | -9.0               | -2.8               |

The table illustrates that the gas concentrations in the alveoli, where the breathing gas is saturated with water vapor at body temperature (BTPS), are 6.2% lower than the corresponding concentrations in the same gas mixture after removal of all water vapor (ATPD).

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| <b>Multi-gas, Sidestream (MASIMO SWEDEN AB)</b>                             |                                                                                                     |
|-----------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------|
| ISA CO2                                                                     | CO2, CO2 waveform                                                                                   |
| ISA AX+                                                                     | CO2, N2O, primary and secondary Agents (HAL, ISO, ENF, SEV, DES)                                    |
| ISA OR+                                                                     | CO2,O2, N2O, primary and secondary Agents (HAL, ISO, ENF, SEV, DES)                                 |
| Gas /CO2 Interface                                                          | Connector and S/W Interface Driver, Applicable for All Gas and CO2 Modules.                         |
| Description                                                                 | Compact, low-flow sidestream gas analyzers with integrated pump, zeroing valve and flow controller. |
| Ambient CO2                                                                 | $\leq 800$ ppm (0.08 vol%)                                                                          |
| Recovery time after defibrillator test                                      | Unaffected                                                                                          |
| Water handling                                                              | Nomoline Family sampling lines with proprietary water removal tubing.                               |
| Sampling flow rate                                                          | $50 \pm 10$ sml/min <sup>1</sup>                                                                    |
| Degree of protection against harmful ingress of water or particulate matter | IPX4                                                                                                |
| Method of sterilization                                                     | The ISA system contains no sterile parts.                                                           |
| Mode of operation                                                           | CONTINUOUS OPERATION                                                                                |
| Degree of protection against electric shock                                 | Nomoline Family sampling lines are classified as DEFIBRILLATION-PROOF TYPE BF APPLIED PART          |
| <b>Data output</b>                                                          |                                                                                                     |
| Breath detection                                                            | Adaptive threshold, minimum 1 vol% change in CO2concentration.                                      |
| Respiration rate <sup>2</sup>                                               | 0 to $150 \pm 1$ breaths/min (or BrPM)                                                              |

<sup>1</sup> Volumetric flow rate of air corrected to standardized conditions of temperature and pressure.  
<sup>2</sup> Measured at I/E ratio 1:1 using breath simulator according to EN ISO 80601-2-55 fig. 201.101.

## Fi and ET<sup>1</sup>

Fi and ET are displayed after one breath and have a continually updated breath average.

The following methods are used to calculate end-tidal (ET) values:

-CO<sub>2</sub>: The highest concentration of CO<sub>2</sub> during one breathing cycle with a weight function applied to favor values closer to the end of the cycle.

-O<sub>2</sub>: The highest/lowest concentration of O<sub>2</sub> during the expiratory phase (depending on whether ET<sub>O2</sub> is higher or lower than Fi<sub>O2</sub>)

-N<sub>2</sub>O and anesthetic agents: The momentary gas concentration at the time point where ET<sub>CO2</sub> is detected.

ET will typically decrease below nominal value (ET<sub>nom</sub>) when respiration rate (RR) exceeds the RR threshold (RR<sub>th</sub>) according to the following formulas:

|                                |                                                                                                           |                                                                                                   |
|--------------------------------|-----------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------|
| ISA CO2                        | ET=ETnom×(125/RR)_for RRth >125                                                                           |                                                                                                   |
| CO2                            |                                                                                                           |                                                                                                   |
| ISA OR+/AX+                    |                                                                                                           |                                                                                                   |
| CO2                            | ET=ETnom×√(70 / RR) for RRth >70                                                                          |                                                                                                   |
| N2O, O2, DES, ENF, ISO, SEV    | ET=ETnom×√(50 / RR) for RRth >50                                                                          |                                                                                                   |
| HAL                            | ET=ETnom×√(35 / RR) for RRth >35                                                                          |                                                                                                   |
| Automatic agent identification | ISA OR+/AX+: primary and secondary agent.                                                                 |                                                                                                   |
| Gas analyzer                   |                                                                                                           |                                                                                                   |
| Sensor head                    | 2 to 9 channel NDIR type gas analyzer measuring at 4 to 10 μm.                                            |                                                                                                   |
|                                | Data acquisition rate 10 kHz (sample rate 20 Hz / channel).                                               |                                                                                                   |
|                                | O2 measurements by Servomex's paramagnetic sensor.                                                        |                                                                                                   |
| Calibration                    | No span calibration is required for the IR bench. An automatic zeroing is performed 1 to 3 times per day. |                                                                                                   |
| Compensation                   | ISA CO2                                                                                                   | Automatic compensation for pressure and temperature.                                              |
|                                |                                                                                                           | Manual compensation for broadening effects on CO2.                                                |
|                                | ISA OR+/AX+                                                                                               | Automatic compensation for pressure, temperature and broadening effects on CO2.                   |
| Warm-up time                   | ISA CO2:                                                                                                  | <10 seconds (concentrations reported and full accuracy)                                           |
|                                | ISA OR+/AX                                                                                                | +: <20 seconds (concentrations reported, automatic agent identification enabled and full accuracy |

Measured at I/E ratio 1:1 using breath simulator according to EN ISO 80601-2-55 fig. 201.101.<sup>1</sup>

|                                                                                                                    |                                                                                              |                                                                 |
|--------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------|-----------------------------------------------------------------|
| Rise time <sup>1</sup>                                                                                             | CO2                                                                                          | ≤ 200 ms(≤ 300 ms for ISA OR+/AX+)                              |
|                                                                                                                    | N2O, O2, ENF, ISO, SEV, DES                                                                  | ≤ 400 ms                                                        |
|                                                                                                                    | HAL                                                                                          | ≤ 500 ms                                                        |
| Primary agent threshold (ISA OR+/AX+)                                                                              | 0.15 vol%. When an agent is identified, concentrations will be reported even below 0.15 vol% |                                                                 |
| Secondary agent threshold (ISA OR+/AX+)                                                                            | 0.2 vol% +10% of total agent concentration                                                   |                                                                 |
| Agent identification time (ISA OR+/AX+)                                                                            | <20 seconds (typically <10 seconds)                                                          |                                                                 |
| Total system response time <sup>2</sup>                                                                            | ISA CO2:                                                                                     | < 3 seconds                                                     |
|                                                                                                                    | ISA OR+/AX+:                                                                                 | < 4 seconds (with 2m Nomoline Airway Adapter Set sampling line) |
| Accuracy standard conditions                                                                                       |                                                                                              |                                                                 |
| The following accuracy specifications are valid with no drift for dry single gases at 22 ± 5 °C and 1013 ± 40 hPa: |                                                                                              |                                                                 |
|                                                                                                                    | Range <sup>3</sup>                                                                           | Accuracy                                                        |
| CO2                                                                                                                | 0 to 15 vol%                                                                                 | ±(0.2 vol% +2% of reading)                                      |
|                                                                                                                    |                                                                                              |                                                                 |
| N2O                                                                                                                | 0 to 100 vol%                                                                                | ±(2 vol% +2% of reading)                                        |
| HAL, ENF, ISO                                                                                                      |                                                                                              | ±(0.15 vol%+5% of reading)                                      |
|                                                                                                                    | 8 to 25 vol%                                                                                 |                                                                 |
| SEV                                                                                                                | 0 to 10 vol%                                                                                 | ±(0.15 vol% +5% of reading)                                     |
|                                                                                                                    |                                                                                              |                                                                 |
| DES                                                                                                                | 0 to 22 vol%                                                                                 | ±(0.15 vol% +5% of reading)                                     |
|                                                                                                                    |                                                                                              |                                                                 |
| O2                                                                                                                 | 0 to 100 vol%                                                                                | ±(1 vol% +2% of reading)                                        |
| Accuracy - all conditions                                                                                          |                                                                                              |                                                                 |

Measured according to EN ISO 80601-2-55.<sup>1</sup>

Measured according to EN ISO 80601-2-55.<sup>2</sup>

All gas concentrations are reported in units of volume percent and may be translated into mmHg or kPa by using the reported atmospheric<sup>3</sup> pressure.

The following accuracy specifications are valid with no drift for all specified environmental conditions, except for interference from water vapor in the below section “Effects from water vapor partial pressure on gas readings”.

| GAS                 | Accuracy                                         |
|---------------------|--------------------------------------------------|
| CO2                 | $\pm(0.3 \text{ kPa} + 4\% \text{ of reading})$  |
| N2O                 | $\pm(2 \text{ kPa} + 5\% \text{ of reading})$    |
| Agents <sup>1</sup> | $\pm(0.2 \text{ kPa} + 10\% \text{ of reading})$ |
| O2                  | $\pm(2 \text{ kPa} + 2\% \text{ of reading})$    |

#### Effects from water vapor partial pressure on gas readings

When the breathing gas flows through the sampling line, the gas temperature will adapt to the ambient temperature before reaching the gas analyzer. The measurement of all gases will always show the actual partial pressure at the current humidity level in the gas sample. As the NOMO section removes all condensed water, no water will reach the ISA gas analyzer. However at an ambient temperature of 37 °C and a breathing gas with a relative humidity of 95% the gas reading will typically be 6% lower than corresponding partial pressure after removal of all water.

<sup>1</sup> The accuracy specification is not valid if more than two agents are present in the gas mixture. If more than two agents are present, an alarm will be set.

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| <b>BFA (Brain Function Assessment)</b> |                                                              |            |
|----------------------------------------|--------------------------------------------------------------|------------|
| BFA Interface                          | Required for Integrating BFA module and monitors             |            |
| EEG sensitivity                        | $\pm 450\mu V$                                               |            |
| Noise                                  | $< 2\mu V_{p-p}$ $< 0.4\mu V$ RMS, 0.25-250 Hz               |            |
| CMRR                                   | $> 140dB$                                                    |            |
| Input impedance                        | $> 50M\Omega$                                                |            |
| Sample rate                            | 1000 samples/sec(16 bits equivalent)                         |            |
| Brain Function Index (BFI)             | 0-100. Filter 1-47Hz, 1sec. update                           |            |
| EMG                                    | 0-100. Filter 30-47 Hz, 1 sec. update                        |            |
| BSR                                    | 0-100. Filter 2-47 Hz, 1 sec. update                         |            |
| Signal Quality Index (SQI)             | 0-100. 1 sec. update                                         |            |
| EEG Waveform                           | $\pm 250\mu V$ , user-adjustable, 5 sec                      |            |
| Alarms                                 | Auditory and visual, user-adjustable limits                  |            |
| Artifact rejection                     | Automatic                                                    |            |
| Sensor impedance measurement           | 0-30kOhm / Manual-Automatic/ measurement current $0.06\mu A$ |            |
| Power supply                           | 5 VDC                                                        |            |
| Power Consumption                      | Less than 0.5 W                                              |            |
| Weight                                 | 100 gr                                                       |            |
| Dimensions                             | 111×64×25 mm                                                 |            |
| Classification                         | Class I, type BF, continuous use                             |            |
| Sensors                                | Ambu Neuro Sensors                                           |            |
| Cable length                           | 195 cm/ 77" with 35 cm/ 14" split                            |            |
| Memory                                 | Data recording (96 hours)                                    |            |
| Trend                                  | BFI/EMG/SQI/BS, 10 sec. update                               |            |
| Environment - Operation                | Temperature                                                  | 5-40°C     |
|                                        | Rel humidity                                                 | 20~96%     |
|                                        | Altitude                                                     | -200~3000m |

| <b>Recorder</b> |                                                          |
|-----------------|----------------------------------------------------------|
| Model           | SAADAT Thermal Printer                                   |
| Channel         | Aria:<br>Up to 2 waveforms<br><br>Aria-TC:<br>1 waveform |
| Printing Speed  | 6, 12.5, 25 mm/sec                                       |
| Paper Size      | 57mm                                                     |

| <b>ALARM</b> |                                                                      |
|--------------|----------------------------------------------------------------------|
| Sources      | Error messages, All other parameter limits                           |
| Alarm On/Off | Selectable for all parameters                                        |
| Alert        | Blinking on Display, Volume Selectable Audio Alarms, Light indicator |

| <b>TREND</b>        |                                                                                                                                                                                                                                                 |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Sources             | Aria:<br>HR,ST,PVCs,AFIB, RESP, T1,T2, IBP1(SYS,DIA,MAP), IBP2(SYS,DIA,MAP)<br>IBP3(SYS,DIA,MAP)<br>IBP4(SYS,DIA,MAP)<br>SPO2, PR,SpHb, PI, SpCo, SpMet, PVI, SpOc,<br>EtCo2,FiCo2,AWRR<br><br>Aria-TC:<br>HR, SpO2, PR, RR, T1, PVCs, ST, AFIB |
| Trend Time Save     | 96 Hours                                                                                                                                                                                                                                        |
| Trend Time Interval | 5, 10, 15, 30, 45 Min, 1, 2, 4 Hours                                                                                                                                                                                                            |
| Resolution          | 1 sec                                                                                                                                                                                                                                           |

| <b>INPUT/OUTPUT</b> |                                                                                                                                                                                            |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Network             | Aria: Digital, TCP/IP (Wi-Fi) and TCP/IP (Wire)<br><br>Aria-TC: TCP/IP (wifi 2.4 GHz), 3g/4g modem (GPRS: HSPA/UMTS: 1900/2100MHz<br>GSM/GPRS/EDGE: 850/900/1800/1900MHz), GSM 0.9/1.8 GHz |

| Internal Battery                        |                   |                      |                                                                                                                                |
|-----------------------------------------|-------------------|----------------------|--------------------------------------------------------------------------------------------------------------------------------|
| Nickel-Metal Hybride 3.6V,2.5AH         |                   |                      |                                                                                                                                |
| Lithium Polymer 11.1V,4.3AH (Aria only) |                   |                      |                                                                                                                                |
| Lithium ion 11.1V,3.3AH                 |                   |                      |                                                                                                                                |
| System Model                            |                   | Nickel-Metal Hybride |                                                                                                                                |
|                                         |                   | Charge time          | Usage                                                                                                                          |
| ARIA                                    |                   | ~ 3hours             | ~ 2:30hours                                                                                                                    |
| F1                                      |                   |                      |                                                                                                                                |
| F1R                                     |                   |                      |                                                                                                                                |
| System Model                            |                   | Lithium Polymer      |                                                                                                                                |
|                                         |                   | Charge time          | Usage                                                                                                                          |
| ARIA                                    |                   |                      |                                                                                                                                |
| F1                                      |                   | ~ 6hours             | ~ 5hours                                                                                                                       |
| F1R                                     |                   | ~ 6hours             | ~ 4hours                                                                                                                       |
| System Model                            |                   | Lithium ion          |                                                                                                                                |
|                                         |                   | Charge time          | Usage                                                                                                                          |
| ARIA                                    |                   |                      |                                                                                                                                |
| F1                                      |                   | ~ 6hours             | ~ 10hours                                                                                                                      |
| F1R                                     |                   | ~ 6hours             | ~ 8hours                                                                                                                       |
| TC station                              | <b>Li-ION2200</b> | ~ 4:30hours          | Normal operation= ~ 3hours<br>In full Load(NIBP=30,12 Leads Recorder & Data Send=25,BackLight=Max,SPO2 Connected)= ~ 1:30hours |
|                                         | <b>Li-ION2600</b> | ~ 6hours             | Normal operation= ~ 4hours<br>In full Load(NIBP=45,12 Leads Recorder & Data Send=30,BackLight=Max,SPO2 Connected)= ~ 2:30hours |
|                                         | <b>Li-ION3200</b> | ~ 6hours             | Normal operation= ~ 5hours<br>In full Load(NIBP=70,12 Leads Recorder & Data Send=50,BackLight=Max,SPO2 Connected)= ~ 3hours    |

| Physical Specification |                         |                        |
|------------------------|-------------------------|------------------------|
| Model                  | Dimension (mm)          | Weight (approximately) |
| ARIA Monitor           | 155(W) × 107(H) × 65(D) | Less than 800g         |
| F1 Station             | 190(W) × 155(H) × 80(D) | 800g                   |
| F1R Station            | 220(W) × 155(H) × 90(D) | 1100g                  |
| ARIA With TC Station   | 235(W) × 225(H) × 90(D) | Less than 3Kg          |

| ENVIRONMENTAL |                                      |                                        |
|---------------|--------------------------------------|----------------------------------------|
| Temperature   | Operating: °C                        | Aria: 5 to 40°C<br>Aria-TC: 0 to 55 °C |
|               | Storage & Transport:                 | -25 to 60 °C                           |
| Humidity      | Operating: (Noncondensing)           | 20-90 % (Noncondensing)                |
|               | Storage & Transport: (Noncondensing) | 10-100 %                               |
| Altitude      | -200 to 3000 m                       |                                        |

| Viewer Specification                     |                                                   |
|------------------------------------------|---------------------------------------------------|
| Storage                                  |                                                   |
| 12 Lead ECG Signal                       | 1000 Records                                      |
| Physician Measurement and Interpretation | 1000 Records                                      |
| Physiological Parameters                 | 1000 Records                                      |
| Print                                    |                                                   |
| Laser Printer                            | Print in Any Size paper                           |
| File                                     | PDF/ JPEG Format                                  |
| Filters                                  |                                                   |
| Notch Filter                             | 50/60 Hz                                          |
| Drift Filter                             | 0.5 Hz                                            |
| Display                                  |                                                   |
| 12 Lead ECG Signal                       | I, II, III, aVR, aVL, aVF, V1, V2, V3, V4, V5, V6 |
| ECG Sample Rate                          | 500                                               |
| ECG Symbol Length                        | 16 bit                                            |
| ECG Signal Length                        | 10 Sec                                            |
| Physiological Parameters                 | Heart Rate, NIBP, SPO2, TEMP1, RR                 |
| Calibration Signal                       | 1 mV, 200 ms                                      |

|                               |                                         |
|-------------------------------|-----------------------------------------|
| Manually Lead selection       | Yes                                     |
| Superimposition               | Yes                                     |
| Patient information           | Name, Patient ID, Gender, Age           |
| Sender Information            | Ambulance ID                            |
| Time Sweep                    | (12.5/25/50) mm/Sec                     |
| Voltage Gain                  | (5/10/20/40) mm/mV                      |
| Physician information         | Name, ID, Interpretation Note           |
| <b>Measurement</b>            |                                         |
| Automatic Measurement         | Optional                                |
| Manual Measurement            | P and QRS Duration, PQ and QT Intervals |
| Heart Axis                    | P, QRS, T Axis                          |
| <b>Other</b>                  |                                         |
| Portable Software             | Yes                                     |
| Touch Screen                  | Yes                                     |
| Compatibility                 | Win XP/Vista/7/8/10                     |
| Upgrade Capability            | Manual                                  |
| <b>Connection</b>             |                                         |
| Connecting to Data Repository | Online/ Offline                         |
| Data Base                     | Costume Format                          |

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| Menu item                          | selection                                                                | Default                             |
|------------------------------------|--------------------------------------------------------------------------|-------------------------------------|
| <b>The parameters in ECG menu</b>  |                                                                          |                                     |
| ECG LEAD                           | I,II,III,aVR,aVF,aVL,V1,V2,V3,V4,V5,V6                                   | II                                  |
| CABLE TYPE                         | 3 Wires,5 Wires, 10 Wires                                                | 3 Wires (10 Wires for Aria-TC)      |
| DISPLAY FORMAT                     | Cascade, 2Traces, 4Traces, 7Traces, 12Traces                             | Cascade                             |
| ECG SIZE                           | CHANGE ( $\times 0.25, \times 0.5, \times 1, \times 2, \times 4$ ), AUTO | AUTO                                |
| ECG SWEEP                          | 12.5, 25, 50mm/s                                                         | 25                                  |
| ECG FILTER                         | MONITOR, NORMAL, EXTENDED                                                | NORMAL                              |
| HR AVERAGE                         | 4, 8, 16SEC                                                              | 8SEC                                |
| HR SOURCE                          | ECG, SpO2, IBP1, IBP2, AUTO                                              | AUTO                                |
| BEAT VOLUME                        | 1, 2, 3, 4, 5, 6, 7, OFF                                                 | 1                                   |
| PACE DETECT                        | ON, OFF                                                                  | OFF                                 |
| ECG CALIB                          | ON, OFF                                                                  | OFF                                 |
| ALARM LEVEL                        | 1, 2                                                                     | 1                                   |
| HR ALARM                           | ON, OFF                                                                  | OFF                                 |
| HR HIGH ALARM                      | HR LOW ALARM +5 to 250                                                   | 150Bpm                              |
| HR LOW ALARM                       | 30 to HR HIGH ALARM -5                                                   | 50Bpm                               |
| <b>The parameters in RESP menu</b> |                                                                          |                                     |
| Menu item                          | selection                                                                | Default                             |
| RESP LEAD.                         | RA-LA, RA-LL                                                             | RA-LA                               |
| RESP GAIN                          | $\times 0.25, \times 0.5, \times 1, \times 2, \times 4$                  | $\times 1$                          |
| RESP SWEEP                         | 3, 6, 12.5, 25mm/s                                                       | 6mm/s                               |
| ALARM LEVEL                        | 1, 2                                                                     | 1                                   |
| RR ALARM                           | ON, OFF                                                                  | OFF                                 |
| RR HIGH ALARM                      | RR LOW ALARM +1 to 150                                                   | 25Brpm                              |
| RR LOW ALARM                       | 5 to RR HIGH ALARM -1                                                    | 5Brpm                               |
| APNEA LIMIT                        | 10 to 40S, OFF                                                           | 10S                                 |
| <b>The parameters in SpO2 menu</b> |                                                                          |                                     |
| Menu item                          | selection                                                                | Default                             |
| AVERAGE TIME                       | 2-4, 4-6, 8, 10, 12, 14, 16                                              | 8                                   |
| SpO2 PLETH SWEEP                   | 12.5, 25mm/s                                                             | Aria: 25 mm/s<br>Aria-TC: 12.5 mm/s |
| ALARM LEVEL                        | 1, 2                                                                     | 1                                   |
| SpO2 ALARM                         | ON, OFF                                                                  | OFF                                 |
| SpO2 HIGH ALARM                    | SpO2 LOW ALARM +1 to 100 (with step 1)                                   | 100                                 |
| SpO2 LOW ALARM                     | 1 to SpO2 HIGH ALARM -1 (with step 1)                                    | 90                                  |
| PR HIGH ALARM                      | PR LOW ALARM +5 to 235                                                   | 140                                 |
| PR LOW ALARM                       | 20 to PR HIGH ALARM -5                                                   | 50                                  |
| SpO2 SENSITVITY MODE               | NORMAL, MAX, APOD                                                        | NORMAL                              |
| SPO2 PULSE RATE                    | ON, OFF                                                                  | ON                                  |
| PI HIGH ALARM                      | PI LOW ALARM +0.1 to 19.0 (with step 0.1)                                | 19.0                                |
| PI LOW ALARM                       | 0.0 to PI HIGH ALARM -0.1 (with step 0.1)                                | 0.0                                 |
| PVI HIGH ALARM                     | PVI LOW ALARM +1 to 99 (with step 1)                                     | 99                                  |
| PVI LOW ALARM                      | 1 to PVI HIGH ALARM -1 (with step 1)                                     | 1                                   |
| SpOC HIGH ALARM                    | SpOC LOW ALARM +1 to 34.0 (with step 1)                                  | 34.0                                |
| SpOC LOW ALARM                     | 1.0 to SpOC HIGH ALARM -1 (with step 1)                                  | 1.0                                 |

|                                    |                                                                                                                              |                                                         |
|------------------------------------|------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------|
| <b>SpCO HIGH ALARM</b>             | SpCO LOW ALARM +1 to 99.0 (with step 1)                                                                                      | 10.0                                                    |
| <b>SpCO LOW ALARM</b>              | 1.0 to SpCO HIGH ALARM -1 (with step 1)                                                                                      | 1.0                                                     |
| <b>SpMet HIGH ALARM</b>            | SpMet LOW ALARM +0.5 to 99.5 (with step 0.5)                                                                                 | 3.0                                                     |
| <b>SpMet LOW ALARM</b>             | 0.5 to SpMet HIGH ALARM -0.5 (with step 0.5)                                                                                 | 0.5                                                     |
| <b>SpHb HIGH ALARM</b>             | SpHb LOW ALARM +0.1 to 24.5 (with step 0.1)                                                                                  | 17.0                                                    |
| <b>SpHb LOW ALARM</b>              | 0.5 to SpHb HIGH ALARM -0.1 (with step 0.1)                                                                                  | 7.0                                                     |
| <b>The parameters in NIBP menu</b> |                                                                                                                              |                                                         |
| <b>Menu item</b>                   | <b>selection</b>                                                                                                             | <b>Default</b>                                          |
| <b>NIBP UNIT</b>                   | mmHg , KPa                                                                                                                   | mmHg                                                    |
| <b>ALARM LEVEL</b>                 | 1,2                                                                                                                          | 1                                                       |
| <b>NIBP ALARM</b>                  | ON,OFF                                                                                                                       | OFF                                                     |
| <b>SYS HIGH ALARM</b>              | Adult: SYS LOW ALARM +5 to 255<br>Neonate: SYS LOW ALARM +5 to 135<br>Pediatric: SYS LOW ALARM +5 to 240<br>(with step 5)    | Adult: 160mmHg<br>Neonate: 90mmHg<br>Pediatric: 120mmHg |
| <b>SYS LOW ALARM</b>               | Adult: 30 to SYS HIGH ALARM -5<br>Neonate: 30 to SYS HIGH ALARM -5<br>Pediatric: 30 to SYS HIGH ALARM -5<br>(with step 5)    | Adult: 90mmHg<br>Neonate: 40mmHg<br>Pediatric: 70mmHg   |
| <b>DIA HIGH ALARM</b>              | Adult: DIA LOW ALARM +5 to 220<br>Neonate: DIA LOW ALARM +5 to 110<br>Pediatric: DIA LOW ALARM +5 to 220<br>(with step 5)    | Adult: 90mmHg<br>Neonate: 60mmHg<br>Pediatric: 70mmHg   |
| <b>DIA LOW ALARM</b>               | Adult: 15 to DIA HIGH ALARM -5<br>Neonate: 15 to DIA HIGH ALARM -5<br>Pediatric: 15 to DIA HIGH ALARM -5<br>(with step 5)    | Adult: 50mmHg<br>Neonate: 20mmHg<br>Pediatric: 40mmHg   |
| <b>MAP HIGH ALARM</b>              | Adult: MAP LOW ALARM +5 to 235<br>Neonate: MAP LOW ALARM +5 to 125<br>Pediatric: MAP LOW ALARM +5 to 230<br>(with step 5)    | Adult: 110mmHg<br>Neonate: 70mmHg<br>Pediatric: 90mmHg  |
| <b>MAP LOW ALARM</b>               | Adult: 20 to MAP HIGH ALARM -5<br>Neonate: 20 to MAP HIGH ALARM -5<br>Pediatric: 20 to MAP HIGH ALARM -5<br>(with step 5)    | Adult: 60mmHg<br>Neonate: 25mmHg<br>Pediatric: 50mmHg   |
| <b>AUTO/MANUAL</b>                 | 1min, 2min, 3min, 5min, 10min, 15min, 20min, 30min, 45min, 60min, 90min, 2hr, 4hr, 8hr, 12hr, 16hr, 20hr, 24hr, MANUAL, STAT | MANUAL                                                  |
| <b>The parameters in TEMP menu</b> |                                                                                                                              |                                                         |
| <b>Menu item</b>                   | <b>selection</b>                                                                                                             | <b>Default</b>                                          |
| <b>TEMP UNIT</b>                   | °C, °F                                                                                                                       | °C                                                      |
| <b>ALARM LEVEL</b>                 | 1,2                                                                                                                          | 1                                                       |
| <b>TEMP ALARM</b>                  | ON ,OFF                                                                                                                      | OFF                                                     |
| <b>T1 HIGH ALARM</b>               | T1 LOW ALARM +0.5 to 50                                                                                                      | 39                                                      |
| <b>T1 LOW ALARM</b>                | 0 to T1 HIGH ALARM -1                                                                                                        | 35                                                      |
| <b>T2 HIGH ALARM</b>               | T2 LOW ALARM +0.5 to 50                                                                                                      | 40                                                      |
| <b>T2 LOW ALARM</b>                | 0 to T2 HIGH ALARM -1                                                                                                        | 36                                                      |
| <b>DT HIGH ALARM</b>               | DT LOW ALARM +1 to 50                                                                                                        | 5                                                       |
| <b>DT LOW ALARM</b>                | 0 to DT HIGH ALARM -1                                                                                                        | 1.0                                                     |
| <b>The parameters in IBP menu</b>  |                                                                                                                              |                                                         |
| <b>Menu item</b>                   | <b>selection</b>                                                                                                             | <b>Default</b>                                          |

|                                     |                                            |                                                  |
|-------------------------------------|--------------------------------------------|--------------------------------------------------|
| <b>IBP UNIT</b>                     | mmHg , KPa,cmH2O                           | mmHg                                             |
| <b>IBP LABEL</b>                    | IBP, ART, PAP, CVP, LAP, RAP, LVP, RVP,ICP | IBP                                              |
| <b>IBP SWEEP</b>                    | 3,6,12.5,25 mm/s                           | 12.5 mm/s                                        |
| <b>IBP GRID</b>                     | ON, OFF                                    | OFF                                              |
| <b>IBP FILTER</b>                   | 8, 16, 22 Hz                               | 16 Hz                                            |
| <b>ALWAYS AUTO SCALE</b>            | ON,OFF                                     | OFF                                              |
| <b>IBP ALARM</b>                    | ON,OFF                                     | OFF                                              |
| <b>ART CATH.<br/>DISCONNECT ALM</b> | ON , OFF                                   | OFF                                              |
| <b>PPV MEASUREMENT</b>              | ON , OFF                                   | OFF                                              |
| <b>ALARM LEVEL</b>                  | 1,2                                        | 1                                                |
| <b>IBP HIGH ALARM</b>               | IBP LOW ALARM +5 to 300                    | SYS: 150 mmHg<br>DIA: 100 mmHg<br>MEAN: 115 mmHg |
| <b>IBP LOW ALARM</b>                | -50to IBP HIGH ALARM -5                    | SYS: 80 mmHg<br>DIA: 50 mmHg<br>MEAN: 60 mmHg    |
| <b>ART HIGH ALARM</b>               | ART LOW ALARM +5 to 300                    | SYS: 150 mmHg<br>DIA: 100 mmHg<br>MEAN: 115 mmHg |
| <b>ART LOW ALARM</b>                | -50to ART HIGH ALARM -5                    | SYS: 80 mmHg<br>DIA: 50 mmHg<br>MEAN: 60 mmHg    |
| <b>LVP HIGH ALARM</b>               | LVP LOW ALARM +5 to 300                    | SYS: 150 mmHg<br>DIA: 20 mmHg<br>MEAN: 80 mmHg   |
| <b>LVP LOW ALARM</b>                | -50 to LVP HIGH ALARM -5                   | SYS: 80 mmHg<br>DIA: -5 mmHg<br>MEAN: 20 mmHg    |
| <b>PAP HIGH ALARM</b>               | PAP LOW ALARM +1 to 120                    | SYS: 40 mmHg<br>DIA: 20 mmHg<br>MEAN: 30 mmHg    |
| <b>PAP LOW ALARM</b>                | -50 to PAP HIGH ALARM -1                   | SYS: 5 mmHg<br>DIA: -5 mmHg<br>MEAN: 0 mmHg      |
| <b>RVP HIGH ALARM</b>               | RVP LOW ALARM +1 to 100                    | SYS: 40 mmHg<br>DIA: 15 mmHg<br>MEAN: 30 mmHg    |
| <b>RVP LOW ALARM</b>                | -50 to RVP HIGH ALARM -1                   | SYS: 5mmHg<br>DIA: -5 mmHg<br>MEAN: 0 mmHg       |
| <b>CVP HIGH ALARM</b>               | CVP LOW ALARM +1 to 100                    | 15 mmHg                                          |
| <b>CVP LOW ALARM</b>                | -50 to CVP HIGH ALARM -1                   | -5 mmHg                                          |
| <b>LAP HIGH ALARM</b>               | LAP LOW ALARM +1 to 100                    | 20 mmHg                                          |
| <b>LAP LOW ALARM</b>                | -50 to LAP HIGH ALARM -1                   | -5 mmHg                                          |
| <b>RAP HIGH ALARM</b>               | RAP LOW ALARM +1 to 100                    | 15 mmHg                                          |
| <b>RAP LOW ALARM</b>                | -50 to RAP HIGH ALARM -1                   | -5 mmHg                                          |

|                       |             |                               |                                                     |
|-----------------------|-------------|-------------------------------|-----------------------------------------------------|
|                       |             |                               |                                                     |
| <b>ICP HIGH ALARM</b> |             | ICP LOW ALARM +1 to 100       | Adult: 10mmHg<br>Neonate: 4mmHg<br>Pediatric: 4mmHg |
| <b>ICP LOW ALARM</b>  |             | -40 to ICP HIGH ALARM -1      | 0 mmHg                                              |
| <b>IBP SCALE</b>      |             |                               |                                                     |
| <b>IBP</b>            | <b>HIGH</b> | LOW +10 TO 300 (with step 10) | 200                                                 |
|                       | <b>LOW</b>  | -50 TO HIGH-10                | -20                                                 |
|                       | <b>SIGN</b> | (HIGH+LOW)/2                  | 90                                                  |
| <b>ART</b>            | <b>HIGH</b> | LOW +10 TO 300 (with step 10) | 200                                                 |
|                       | <b>LOW</b>  | -50 TO HIGH-10                | 40                                                  |
|                       | <b>SIGN</b> | (HIGH+LOW)/2                  | 120                                                 |
| <b>PAP</b>            | <b>HIGH</b> | LOW +5 TO 300 (with step 5)   | 80                                                  |
|                       | <b>LOW</b>  | -50 TO HIGH-5                 | -10                                                 |
|                       | <b>SIGN</b> | (HIGH+LOW)/2                  | 35                                                  |
| <b>CVP</b>            | <b>HIGH</b> | LOW +5 TO 300 (with step 5)   | 30                                                  |
|                       | <b>LOW</b>  | -50 TO HIGH-5                 | -10                                                 |
|                       | <b>SIGN</b> | (HIGH+LOW)/2                  | 10                                                  |
| <b>LAP</b>            | <b>HIGH</b> | LOW +5 TO 300 (with step 5)   | 40                                                  |
|                       | <b>LOW</b>  | -50 TO HIGH-5                 | -10                                                 |
|                       | <b>SIGN</b> | (HIGH+LOW)/2                  | 15                                                  |
| <b>RAP</b>            | <b>HIGH</b> | LOW +5 TO 300 (with step 5)   | 30                                                  |
|                       | <b>LOW</b>  | -50 TO HIGH-5                 | -10                                                 |
|                       | <b>SIGN</b> | (HIGH+LOW)/2                  | 10                                                  |
| <b>LVP</b>            | <b>HIGH</b> | LOW +10 TO 300 (with step 10) | 200                                                 |
|                       | <b>LOW</b>  | -50 TO HIGH-10                | -20                                                 |
|                       | <b>SIGN</b> | (HIGH+LOW)/2                  | 90                                                  |
| <b>RVP</b>            | <b>HIGH</b> | LOW +5 TO 300 (with step 5)   | 80                                                  |
|                       | <b>LOW</b>  | -50 TO HIGH-5                 | -10                                                 |
|                       | <b>SIGN</b> | (HIGH+LOW)/2                  | 35                                                  |
| <b>ICP</b>            | <b>HIGH</b> | LOW +5 TO 100 (with step 5)   | 40                                                  |
|                       | <b>LOW</b>  | -40 TO HIGH-5                 | -10                                                 |
|                       | <b>SIGN</b> | (HIGH+LOW)/2                  | 15                                                  |

| The parameters in ARR menu |               |                                             |
|----------------------------|---------------|---------------------------------------------|
| Menu item                  | selection     | Default                                     |
| ARR MONITOR                | ON, OFF       | OFF                                         |
| ALARM LEVEL                | ASYSTOLE      | 1                                           |
|                            | VFIB          | 1                                           |
|                            | VTAC          | 1                                           |
|                            | RUN           | 1, 2, OFF                                   |
|                            | AIVR          | 1, 2, OFF                                   |
|                            | COUPLET       | 1, 2, OFF                                   |
|                            | BIGEMINY      | 1, 2, OFF                                   |
|                            | TRIGEMINY     | 1, 2, OFF                                   |
|                            | TACHY         | 1, 2, OFF                                   |
|                            | BRADY         | 1, 2, OFF                                   |
|                            | AFIB          | 1, 2, OFF                                   |
|                            | PAUS          | 1, 2, OFF                                   |
|                            | FREQUENT PVCs | 1, 2, OFF                                   |
| RATE                       | VTAC          | 100 to 200 (with step 10)                   |
|                            | RUN           | VTAC <sub>rate</sub>                        |
|                            | AIVR          | <VTAC <sub>rate</sub> -1                    |
|                            | TACHY         | 100 to 200 (with step 10)                   |
|                            | BRADY         | 30 to 105 (with step 5)                     |
| COUNT                      | VTAC          | 5 to 12 (with step 1)                       |
|                            | RUN           | 3 to VTAC <sub>count</sub> -1 (with step 1) |
|                            | AIVR          | -                                           |
|                            | FREQUENT PVCs | 1 to 15 (with step 5)                       |
| ARCHIVE                    | ASYSTOLE      | STR, STR/REC, OFF, REC                      |
|                            | VFIB          | STR, STR/REC, OFF, REC                      |
|                            | VTAC          | STR, STR/REC, OFF, REC                      |
|                            | RUN           | STR, STR/REC, OFF, REC                      |
|                            | AIVR          | STR, STR/REC, OFF, REC                      |
|                            | COUPLET       | STR, STR/REC, OFF, REC                      |
|                            | BIGEMINY      | STR, STR/REC, OFF, REC                      |
|                            | TRIGEMINY     | STR, STR/REC, OFF, REC                      |
|                            | TACHY         | STR, STR/REC, OFF, REC                      |
|                            | BRADY         | STR, STR/REC, OFF, REC                      |
|                            | AFIB          | STR, STR/REC, OFF, REC                      |
|                            | PAUS          | STR, STR/REC, OFF, REC                      |
|                            | FREQUENT PVCs | -                                           |

| The parameters in ST menu |                          |         |
|---------------------------|--------------------------|---------|
| Menu item                 | Selection                | Default |
| ST ANALYSIS               | ON, OFF                  | OFF     |
| ST ALARM                  | ON, OFF                  | OFF     |
| ALARM LEVEL               | 1, 2                     | 1       |
| ST LOW ALARM              | -2 to ST HIGH ALARM -0.1 | -0.2    |
| ST HIGH ALARM             | ST LOW ALARM +0.1 to 2   | 0.2     |
| EVENT DURATION            | 15S, 30S, 45S, 60S, OFF  | OFF     |

| The Parameters in GAS WINDOW(Mainstream & Sidestream) |                                            |                            |                                            |
|-------------------------------------------------------|--------------------------------------------|----------------------------|--------------------------------------------|
| Menu item                                             | selection                                  |                            |                                            |
| CO2 UNIT                                              | KPa ,%V ,mmHg                              |                            | mmHg                                       |
| SIGNAL SWEEP                                          | 3mm/s, 6mm/s, 12.5mm/s, 25mm/s             |                            | 12.5mm/s                                   |
| SIGNAL SCALE                                          | CO2                                        | 6%,10%,Auto scale          | 10%                                        |
|                                                       | O2/N2O                                     | 0-50%,0-100%, Auto scale   | 100%                                       |
|                                                       | AA                                         | 1,2,3,5,10,20%, Auto scale | 20%                                        |
| WAVEFORM (Mainstream)                                 | CO2, N2O, AA                               |                            | CO2                                        |
| WAVEFORM (Sidestream)                                 | CO2, O2, N2O, AA                           |                            | CO2                                        |
| O2 COMPENSATE                                         | 1-100 vol%, OFF                            |                            | 21% , AUTO                                 |
| N2O COMPENSATE                                        | 0-100 vol% ( ONLY FOR ISA CO2, IRMA2 CO2 ) |                            | 0%                                         |
| GAS UNIT                                              | KPa ,%V                                    |                            | %V                                         |
| AGENT                                                 | ISO,ENF,HAL,DES,SEV                        |                            | HAL                                        |
|                                                       |                                            |                            | AUTO (For IRMA(AX+) & ISA(OR+) & ISA(AX+)) |
| WORK MODE                                             | MEASURE, STANDBY                           |                            | MEASURE                                    |
| GAS/RESP                                              | GAS, RESP                                  |                            | GAS                                        |
| FIIL SIGNAL                                           | ON,OFF                                     |                            | OFF                                        |
| CO2 ALARM                                             | ON,OFF                                     |                            | OFF                                        |
| N2O ALARM                                             | ON,OFF                                     |                            | OFF                                        |
| AA ALARM                                              | ON,OFF                                     |                            | OFF                                        |
| O2 ALARM                                              | ON,OFF                                     |                            | OFF                                        |
| ALARM LEVEL                                           | 1,2                                        |                            | 2                                          |
| APNEA ALARM                                           | 20s,25s,30s,35s,40s,45s,50s, 55s,60s, OFF  |                            | 20s                                        |
|                                                       |                                            | ADULT/PED                  | NEONATE                                    |
| AWRR LOW                                              | 1~(HIGH-1)                                 | 5 BrPM                     | 15 BrPM                                    |
| AWRR HIGH                                             | (LOW+1)~120                                | 30 BrPM                    | 60 BrPM                                    |
| EtCo2 LOW                                             | 0.4~(HIGH-0.1) (%V)                        | 2.6%V                      |                                            |
| EtCo2 HIGH                                            | (LOW+0.1)~13(%V)                           | 6.5%V                      |                                            |
| FiCo2 HIGH                                            | 0.4~ 13(%V)                                | 1.3%V                      |                                            |
| EtO2,FiO2 LOW (sidestream)                            | 18~(HIGH-1) (%V)                           | 50%                        |                                            |
| EtO2,FiO2 HIGH (sidestream)                           | (LOW+1)~105(%V)                            | 100%                       |                                            |
| EtN2O ,FiN2O LOW                                      | 1~(HIGH-1) (%V)                            | 35%                        |                                            |
| EtN2O HIGH                                            | (LOW+1)~100(%V)                            | 75%                        |                                            |
| FiN2O HIGH                                            | (LOW+1)~82(%V)                             | 75%                        |                                            |
| EtDES ,FiDES LOW                                      | 0.1~(HIGH-0.1) (%V)                        | 5%                         |                                            |
| EtDES ,FiDES HIGH                                     | (LOW+0.1)~18(%V)                           | 10%                        |                                            |
| EtISO ,FiISO LOW                                      | 0.1~(HIGH-0.1) (%V)                        | 0.8%                       |                                            |
| EtISO ,FiISO HIGH                                     | (LOW+0.1)~5(%V)                            | 2%                         |                                            |
| EtENF ,FiENF LOW                                      | 0.1~(HIGH-0.1) (%V)                        | 0.5%                       |                                            |
| EtENF ,FiENF HIGH                                     | (LOW+0.1)~5(%V)                            | 1.5%                       |                                            |
| EtSEV ,FiSEV LOW                                      | 0.1~(HIGH-0.1) (%V)                        | 1%                         |                                            |
| EtSEV ,FiSEV HIGH                                     | (LOW+0.1)~8(%V)                            | 3%                         |                                            |
| EtHAL ,FiHAL LOW                                      | 0.1~(HIGH-0.1) (%V)                        | 0.5%                       |                                            |

|                          |                     |      |
|--------------------------|---------------------|------|
| <b>EtHAL ,FiHAL HIGH</b> | (LOW+0.1)~5(%V)     | 1.5% |
| <b>ZERO</b>              | Only for Mainstream |      |

| The Parameters in BFA WINDOW            |                                                                                                                                                                                                      |                  |
|-----------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------|
| Menu item                               | selection                                                                                                                                                                                            | Default          |
| EEG Gain                                | 25uV,50-250uV                                                                                                                                                                                        | 100uV            |
| BFA ALARM                               | ON,OFF                                                                                                                                                                                               | OFF              |
| BFI LOW                                 | 1~(HIGH-1)                                                                                                                                                                                           | 35%              |
| BFI HIGH                                | (LOW+1)~100                                                                                                                                                                                          | 60%              |
| The Parameters in Cardiac Output WINDOW |                                                                                                                                                                                                      |                  |
| Catheter Type                           | 131HF7,139HF75P,Simulator                                                                                                                                                                            | 131HF7           |
| Temp_Scale                              | 1,2,4                                                                                                                                                                                                | 1                |
| SYSTEM DEFUALT                          |                                                                                                                                                                                                      |                  |
| ALARM VOLUME                            | 1,2,3,4,5,6,7                                                                                                                                                                                        | 1                |
| CALENDAR                                | SOLAR, CHRISTIAN                                                                                                                                                                                     | CHRISTIAN        |
| PATIENT CAT.                            | ADUL,NEONATE,PEDIATRIC                                                                                                                                                                               | ADULT            |
| BED NUMBER                              | 1..150                                                                                                                                                                                               | 01               |
| TOUCH SOUND                             | 1, 2, 3, OFF                                                                                                                                                                                         | 1                |
| BACK LIGHT                              | 1 to 8                                                                                                                                                                                               | 18.5" Monitor: 7 |
|                                         | 1 to 6                                                                                                                                                                                               | 12" Monitor: 5   |
|                                         |                                                                                                                                                                                                      | 10" Monitor: 3   |
|                                         |                                                                                                                                                                                                      | 15" Monitor: 2   |
| BED TO BED                              | DURATION :1,2,3,4,5                                                                                                                                                                                  | 1 Min            |
| Module Color                            |                                                                                                                                                                                                      |                  |
| ECG COLOR                               | Green                                                                                                                                                                                                | GREEN            |
| IBP1 COLOR                              | White, Blue, Brown, Green, Red, Yellow, Cyan, Orange, Cream, Magenta, Light Brown, Light Green, Light Yellow, Light Red, Light Blue, Light Cyan, Light Orange, Light Magenta, Dark Orange, Dark Cyan | LIGHT RED        |
| IBP2 COLOR                              | White, Blue, Brown, Green, Red, Yellow, Cyan, Orange, Cream, Magenta, Light Brown, Light Green, Light Yellow, Light Red, Light Blue, Light Cyan, Light Orange, Light Magenta, Dark Orange, Dark Cyan | LIGHT BLUE       |
| IBP3 COLOR                              | White, Blue, Brown, Green, Red, Yellow, Cyan, Orange, Cream, Magenta, Light Brown, Light Green, Light Yellow, Light Red, Light Blue, Light Cyan, Light Orange, Light Magenta, Dark Orange, Dark Cyan | DARK ORANGE      |

|                   |                                                                                                                                                                                                      |           |
|-------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------|
| <b>IBP4 COLOR</b> | White, Blue, Brown, Green, Red, Yellow, Cyan, Orange, Cream, Magenta, Light Brown, Light Green, Light Yellow, Light Red, Light Blue, Light Cyan, Light Orange, Light Magenta, Dark Orange, Dark Cyan | DARK CYAN |
| <b>SpO2 COLOR</b> | White, Blue, Brown, Green, Red, Yellow, Cyan, Orange, Cream, Magenta, Light Brown, Light Green, Light Yellow, Light Red, Light Blue, Light Cyan, Light Orange, Light Magenta, Dark Orange, Dark Cyan | MAGENTA   |
| <b>CO2 COLOR</b>  | White, Blue, Brown, Green, Red, Yellow, Cyan, Orange, Cream, Magenta, Light Brown, Light Green, Light Yellow, Light Red, Light Blue, Light Cyan, Light Orange, Light Magenta, Dark Orange, Dark Cyan | YELLOW    |
| <b>RESP COLOR</b> | White, Blue, Brown, Green, Red, Yellow, Cyan, Orange, Cream, Magenta, Light Brown, Light Green, Light Yellow, Light Red, Light Blue, Light Cyan, Light Orange, Light Magenta, Dark Orange, Dark Cyan | YELLOW    |
| <b>NIBP COLOR</b> | White, Blue, Brown, Green, Red, Yellow, Cyan, Orange, Cream, Magenta, Light Brown, Light Green, Light Yellow, Light Red, Light Blue, Light Cyan, Light Orange, Light Magenta, Dark Orange, Dark Cyan | WHITE     |
| <b>TEMP COLOR</b> | White, Blue, Brown, Green, Red, Yellow, Cyan, Orange, Cream, Magenta, Light Brown, Light Green, Light Yellow, Light Red, Light Blue, Light Cyan, Light Orange, Light Magenta, Dark Orange, Dark Cyan | CYAN      |

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| Guidance and manufacturer's declaration – electromagnetic emissions                                                                                                                                            |            |                                                                                                                                                                                                                              |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The <b>ARIA</b> Patient Care Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the <b>ARIA</b> should assure that it is used in such an environment. |            |                                                                                                                                                                                                                              |
| Emissions test                                                                                                                                                                                                 | Compliance | Electromagnetic environment - guidance                                                                                                                                                                                       |
| RF emissions<br>CISPR 11                                                                                                                                                                                       | Group 1    | The <b>ARIA</b> uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.                                         |
| RF emissions<br>CISPR 11                                                                                                                                                                                       | Class B    | The <b>ARIA</b> is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions<br>IEC 61000-3-2                                                                                                                                                                            | Complies   | The <b>ARIA</b> is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Voltage fluctuations/ flicker emissions<br>IEC 61000-3-3                                                                                                                                                       | Complies   |                                                                                                                                                                                                                              |

| Guidance and manufacturer's declaration – electromagnetic immunity                                                                                                                                             |                                                              |                  |                                                                                                                                               |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------|------------------|-----------------------------------------------------------------------------------------------------------------------------------------------|
| The <b>ARIA</b> Patient Care Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the <b>ARIA</b> should assure that it is used in such an environment. |                                                              |                  |                                                                                                                                               |
| Immunity test                                                                                                                                                                                                  | IEC 60601 test level                                         | Compliance level | Electromagnetic environment - guidance                                                                                                        |
| Electrostatic discharge (ESD)<br>IEC 61000-4-2                                                                                                                                                                 | ±6 kV contact<br>±8 kV air                                   | Complies         | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. |
| Electrical fast transient/burst<br>IEC 61000-4-4                                                                                                                                                               | ±2 kV for power supply lines<br>±1 kV for input/output lines | Complies         | Mains power quality should be that of a typical commercial or hospital environment.                                                           |
| Surge<br>IEC 61000-4-5                                                                                                                                                                                         | ±1 kV differential mode<br>±2 kV common mode                 | Complies         | Mains power quality should be that of a typical commercial or hospital environment.                                                           |

|                                                                                      |                                                                                                                                                                                                                          |          |                                                                                                                                                                                                                                                       |
|--------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Voltage dips, short interruptions and voltage variations on power supply input lines | <p>&lt;5% UT<br/>(&gt;95% dip in UT)<br/>for 0.5 cycle</p> <p>40% UT<br/>(60% dip in UT)<br/>for 5 cycles</p> <p>70% UT<br/>(30% dip in UT)<br/>for 25 cycles</p> <p>&lt;5% UT<br/>(&gt;95% dip in UT)<br/>for 5 sec</p> | Complies | Mains power quality should be that of a typical commercial or hospital environment. If the user of the <b>ARIA</b> requires continued operation, it is recommended that the <b>ARIA</b> be powered from an uninterruptible power supply or a battery. |
| Power frequency (50/60 Hz) magnetic field<br>IEC 61000-4-8                           | 3 A/m                                                                                                                                                                                                                    | Complies | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.                                                                                                             |
| NOTE $U_T$ is the a.c. mains voltage prior to application of test level.             |                                                                                                                                                                                                                          |          |                                                                                                                                                                                                                                                       |

| Guidance and manufacturer's declaration – electromagnetic immunity                                                                                                                                                                                 |                                                                        |                         |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------|-------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The <b>ARIA</b> Patient Care Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the <b>ARIA</b> should assure that it is used in such an environment.                                     |                                                                        |                         |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Immunity test                                                                                                                                                                                                                                      | IEC 60601 test level                                                   | Compliance level        | Electromagnetic environment – guidance                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <p>Conducted RF<br/>IEC 61000-4-6</p> <p>Radiated RF<br/>IEC 61000-4-3</p>                                                                                                                                                                         | <p>3 Vrms<br/>150 kHz to 80 MHz</p> <p>3 V/m<br/>80 MHz to 2.5 GHz</p> | <p>3 V</p> <p>3 V/m</p> | <p>Portable and mobile RF communications equipment should be used no closer to any part of the <b>ARIA</b>, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance</b></p> $d = 1.17 \sqrt{P}$ $d = 1.17 \sqrt{P}$ <p>80 MHz to 800 MHz</p> $d = 2.33 \sqrt{P}$ <p>800 MHz to 2.5 GHz</p> <p>Where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in meters (m).<br/>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency range.<sup>b</sup><br/>Interference may occur in the vicinity of equipment marked with the following symbol:</p>  |
| <p>NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p> |                                                                        |                         |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |

- a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the ARIA is used exceeds the applicable RF compliance level above, the ARIA should be observed to verify normal operation. If abnormal performance is observed, additional measures may necessary, such as reorienting or relocating the ARIA.
- b Over the frequency range 150kHz to 80 MHz, field strengths should be less than 3 V/m.

| Recommended separation distances between<br>Portable and mobile RF communications equipment and the Vital Sign Monitor                                                                                                                                                                                                                                                                                                                                            |                                                                |                                             |                                           |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------|---------------------------------------------|-------------------------------------------|
| The <b>ARIA</b> Patient Care Monitor is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the <b>ARIA</b> can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the <b>ARIA</b> as recommended below, according to the maximum output power of the communications equipment. |                                                                |                                             |                                           |
| Rated maximum output power of transmitter<br>W                                                                                                                                                                                                                                                                                                                                                                                                                    | Separation distance according to frequency of transmitter<br>m |                                             |                                           |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                   | 150 kHz to 80<br>MHz<br>$d = 1.17 \sqrt{P}$                    | 80 MHz to 800<br>MHz<br>$d = 1.17 \sqrt{P}$ | 800 MHz to 2.5 GHz<br>$d = 2.33 \sqrt{P}$ |
| 0.01                                                                                                                                                                                                                                                                                                                                                                                                                                                              | 0.12                                                           | 0.12                                        | 0.23                                      |
| 0.1                                                                                                                                                                                                                                                                                                                                                                                                                                                               | 0.37                                                           | 0.37                                        | 0.74                                      |
| 1                                                                                                                                                                                                                                                                                                                                                                                                                                                                 | 1.17                                                           | 1.17                                        | 2.33                                      |
| 10                                                                                                                                                                                                                                                                                                                                                                                                                                                                | 3.70                                                           | 3.70                                        | 7.37                                      |
| 100                                                                                                                                                                                                                                                                                                                                                                                                                                                               | 11.7                                                           | 11.7                                        | 23.3                                      |
| For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.                                                                                                                                  |                                                                |                                             |                                           |
| NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.                                                                                                                                                                                                                                                                                                                                                                     |                                                                |                                             |                                           |
| NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.                                                                                                                                                                                                                                                                                                |                                                                |                                             |                                           |

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| Guidance and manufacturer's declaration – electromagnetic emissions                                                                                                                                                  |            |                                                                                                                                                                                                                                 |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The <b>ARIA TC</b> Patient Care Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the <b>ARIA TC</b> should assure that it is used in such an environment. |            |                                                                                                                                                                                                                                 |
| Emissions test                                                                                                                                                                                                       | Compliance | Electromagnetic environment - guidance                                                                                                                                                                                          |
| RF emissions<br>CISPR 11                                                                                                                                                                                             | Group 2    | The <b>ARIA TC</b> must emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected.                                                                                     |
| RF emissions<br>CISPR 11                                                                                                                                                                                             | Class B    | The <b>ARIA TC</b> is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions<br>IEC 61000-3-2                                                                                                                                                                                  | Complies   | The <b>ARIA TC</b> is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Voltage fluctuations/<br>flicker emissions<br>IEC 61000-3-3                                                                                                                                                          | Complies   |                                                                                                                                                                                                                                 |

| Guidance and manufacturer's declaration – electromagnetic immunity                                                                                                                                                   |                                                                                                                                                                                   |                  |                                                                                                                                                                                                                                                             |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The <b>ARIA TC</b> Patient Care Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the <b>ARIA TC</b> should assure that it is used in such an environment. |                                                                                                                                                                                   |                  |                                                                                                                                                                                                                                                             |
| Immunity test                                                                                                                                                                                                        | IEC 60601 test level                                                                                                                                                              | Compliance level | Electromagnetic environment - guidance                                                                                                                                                                                                                      |
| Electrostatic discharge (ESD)<br>IEC 61000-4-2                                                                                                                                                                       | ±6 kV contact<br>±8 kV air                                                                                                                                                        | Complies         | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.                                                                                                               |
| Electrical fast transient/burst<br>IEC 61000-4-4                                                                                                                                                                     | ±2 kV for power supply lines<br>±1 kV for input/output lines                                                                                                                      | Complies         | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                         |
| Surge<br>IEC 61000-4-5                                                                                                                                                                                               | ±1 kV differential mode<br>±2 kV common mode                                                                                                                                      | Complies         | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                         |
| Voltage dips, short interruptions and voltage variations on power supply input lines                                                                                                                                 | <5% UT<br>(>95% dip in UT) for 0.5 cycle<br><br>40% UT<br>(60% dip in UT) for 5 cycles<br><br>70% UT<br>(30% dip in UT) for 25 cycles<br><br><5% UT<br>(>95% dip in UT) for 5 sec | Complies         | Mains power quality should be that of a typical commercial or hospital environment. If the user of the <b>ARIA TC</b> requires continued operation, it is recommended that the <b>ARIA TC</b> be powered from an uninterruptible power supply or a battery. |
| Power frequency (50/60 Hz) magnetic field<br>IEC 61000-4-8                                                                                                                                                           | 3 A/m                                                                                                                                                                             | Complies         | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.                                                                                                                   |
| NOTE $U_T$ is the a.c. mains voltage prior to application of test level.                                                                                                                                             |                                                                                                                                                                                   |                  |                                                                                                                                                                                                                                                             |
|                                                                                                                                                                                                                      |                                                                                                                                                                                   |                  |                                                                                                                                                                                                                                                             |
|                                                                                                                                                                                                                      |                                                                                                                                                                                   |                  |                                                                                                                                                                                                                                                             |

| Guidance and manufacturer's declaration – electromagnetic immunity                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |                             |                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------|------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The <b>ARIA TC</b> Patient Care Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the <b>ARIA TC</b> should assure that it is used in such an environment.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |                             |                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Immunity test                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              | IEC 60601 test level        | Compliance level | Electromagnetic environment – guidance                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Conducted RF<br>IEC 61000-4-6                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              | 3 Vrms<br>150 kHz to 80 MHz | 3 V              | Portable and mobile RF communications equipment should be used no closer to any part of the <b>ARIA TC</b> , including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.<br><b>Recommended separation distance</b><br>$d = 1.17\sqrt{P}$ $d = 1.17\sqrt{P}$ 80 MHz to 800 MHz<br>$d = 2.33\sqrt{P}$ 800 MHz to 2.5 GHz<br>Where <i>P</i> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <i>d</i> is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range. <sup>b</sup><br>Interference may occur in the vicinity of equipment marked with the following symbol:<br> |
| Radiated RF<br>IEC 61000-4-3                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               | 3 V/m<br>80 MHz to 2.5 GHz  | 3 V/m            |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.<br>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |                             |                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the <b>ARIA TC</b> is used exceeds the applicable RF compliance level above, the <b>ARIA TC</b> should be observed to verify normal operation. If abnormal performance is observed, additional measures may necessary, such as reorienting or relocating the <b>ARIA TC</b> . |                             |                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <sup>b</sup> Over the frequency range 150kHz to 80 MHz, field strengths should be less than 3 V/m.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |                             |                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |

| Recommended separation distances between<br>Portable and mobile RF communications equipment and the <i>Vital Sign Monitor</i>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |                                                                       |                                          |                                           |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------|------------------------------------------|-------------------------------------------|
| The <b>ARIA TC</b> Patient Care Monitor is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the <b>ARIA TC</b> can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the <b>ARIA TC</b> as recommended below, according to the maximum output power of the communications equipment.                                                                                                                                                                       |                                                                       |                                          |                                           |
| Rated maximum output power of transmitter<br><b>W</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            | Separation distance according to frequency of transmitter<br><b>m</b> |                                          |                                           |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  | 150 kHz to 80 MHz<br>$d = 1.17 \sqrt{P}$                              | 80 MHz to 800 MHz<br>$d = 1.17 \sqrt{P}$ | 800 MHz to 2.5 GHz<br>$d = 2.33 \sqrt{P}$ |
| 0.01                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             | 0.12                                                                  | 0.12                                     | 0.23                                      |
| 0.1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              | 0.37                                                                  | 0.37                                     | 0.74                                      |
| 1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                | 1.17                                                                  | 1.17                                     | 2.33                                      |
| 10                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               | 3.70                                                                  | 3.70                                     | 7.37                                      |
| 100                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              | 11.7                                                                  | 11.7                                     | 23.3                                      |
| <p>For transmitters rated at a maximum output power not listed above, the recommended separation distance <math>d</math> in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.</p> <p>NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.</p> <p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p> |                                                                       |                                          |                                           |

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