


bp4 Blood pressure monitor Quick Reference


1. Wrap the cuff on the patient's arm properly.

- Wrap the cuff on a bare arm or over thin clothing. Make sure there is no excessive pressure due to a rolled-up sleeve.
- Place the  sign directly over the brachial artery.
- Place the tail of the cuff within the  sign.
- The bottom edge must be 1 to 2 cm above the inner side of the elbow joint.
- The cuff should be wrapped to a tightness that roughly allows two fingers to be inserted under.



2. Place the patient in proper position.

- It is advised for the patient to rest for 5 minutes prior to measurement.
- The patient's arm must rest on a firm surface.
- Keep the arm on which the cuff is wrapped at the same height as the heart.



3. Start measurement.

- Press the **START/STOP** button to turn the device on.
- Press the **START/STOP** button one more time to start measurement.
- The patient must remain motionless and not speak during measurement.
- When the procedure is finished, the device displays heartrate, systolic and diastolic pressures.
- Hold the **START/STOP** button for 3 seconds to turn the device off.
- The device automatically shuts down if not used for 3 minutes.



For more information about the device, please refer to BP4 User Manual.



Pooyandegan Rah Saadat Co.

After Sales Service: +98 21 73098000, +98 912 1977157

CE 2195

www.saadatco.com

D01201-V1